


# PRO-seq and other methods to study nascent transcription

HDYDI, 9<sup>th</sup> June 2020


Henry F. Thomas

# Why would you care?

- Detailed information about position of transcribing polymerase, TSS, splicing...
- Temporal resolution



Haberle&Stark, 2018



Hirabayashi *et al.*, 2019


# Why would you care?

- Detailed information about position of transcribing polymerase, TSS, splicing...
- Temporal resolution
- Transcriptional vs post-transcriptional regulation




# Principles of nascent RNA isolation

Wissink *et al.*, 2019

## a Chromatin-associated RNA enrichment




## b Pol II-associated RNA enrichment





# Principles of nascent RNA isolation

## c Run-on RNA enrichment



## d Metabolic RNA labelling



# Comparison of methods

Table 1 | Strengths and limitations of nascent RNA methods

Method	Advantages	Considerations
caRNA-seq	<ul style="list-style-type: none"> <li>• Can be used to isolate all chromatin-associated RNA species</li> <li>• Can be combined with methods that assay co-transcriptional processes, including RNA methylation and editing</li> </ul>	Also sequences non-nascent RNAs that stably associate with chromatin
Start-seq	<ul style="list-style-type: none"> <li>• Simultaneously identifies initiation and pausing sites</li> <li>• Allows de novo calling of putative enhancers</li> </ul>	Does not report transcription beyond the first ~100 nucleotides
Yeast NET-seq	<ul style="list-style-type: none"> <li>• Is Pol II specific (antibody enrichment)</li> <li>• Identifies Pol II positions at nucleotide resolution genome-wide</li> </ul>	Is limited to cells with epitope-tagged Pol II
mNET-seq	<ul style="list-style-type: none"> <li>• Is Pol II specific (antibody enrichment)</li> <li>• Identifies Pol II positions at nucleotide resolution genome-wide</li> <li>• Can isolate Pol II with different post-translational modifications</li> </ul>	<ul style="list-style-type: none"> <li>• Includes RNAs that are stably associated with Pol II</li> <li>• Does not currently include RNA &lt;30 nucleotides in length</li> <li>• Has detected eRNA transcription from previously called enhancers</li> </ul>
PRO-cap	<ul style="list-style-type: none"> <li>• Identifies transcription initiation sites</li> <li>• Allows de novo calling of putative enhancers</li> </ul>	Does not report transcription beyond the first ~100 nucleotides
PRO-seq	<ul style="list-style-type: none"> <li>• Captures RNAs from transcriptionally competent polymerases</li> <li>• Identifies positions of active transcription at nucleotide resolution genome-wide</li> <li>• Allows de novo calling of putative enhancers</li> </ul>	<ul style="list-style-type: none"> <li>• Does not measure polymerase backtracking</li> <li>• Also captures RNAs being transcribed from Pol I and Pol III</li> </ul>
CoPRO	<ul style="list-style-type: none"> <li>• Simultaneously identifies initiation and pausing sites</li> <li>• Measures RNA capping status</li> </ul>	Does not measure transcription beyond promoter-proximal pause site
SMIT-seq	Measures splicing status during transcription	Limited to species with short introns
TT-seq	<ul style="list-style-type: none"> <li>• Captures RNAs from actively transcribing polymerases</li> <li>• Can be used to determine RNA stability</li> <li>• Identifies transcription termination sites</li> </ul>	<ul style="list-style-type: none"> <li>• Does not detect Pol II pausing</li> <li>• Has detected eRNA transcription from previously called enhancers</li> </ul>
SLAM-seq and TimeLapse-seq	<ul style="list-style-type: none"> <li>• Captures RNAs from actively transcribing polymerases</li> <li>• Can be used to determine RNA stability</li> </ul>	<ul style="list-style-type: none"> <li>• Requires deep sequencing to measure chemical conversion rate</li> <li>• Long labelling times do not capture newly synthesized RNA</li> </ul>
Intron sequential FISH	<ul style="list-style-type: none"> <li>• Detects transcription of thousands of genes in single cells</li> <li>• Contains positional information of transcribed genes in the 3D space of the nucleus</li> </ul>	<ul style="list-style-type: none"> <li>• Does not report chromosomal positions of active Pol II complexes</li> <li>• Does not distinguish different steps of transcription</li> <li>• Requires a library of intron-targeting probes and series of hybridizations</li> </ul>

caRNA-seq, chromatin-associated RNA sequencing; CoPRO, coordinated precision run-on and sequencing; eRNA enhancer RNA; FISH, fluorescence in situ hybridization; mNET-seq, mammalian native elongating transcript sequencing; NET-seq, native elongating transcript sequencing; Pol, RNA polymerase; PRO-cap, precision run-on with cap selection; PRO-seq, precision run-on sequencing; SLAM-seq, thiol (SH)-linked alkylation for the metabolic sequencing of RNA; SMIT-seq, single-molecule intron tracking sequencing; TT-seq, transient transcriptome sequencing.

Wissink *et al.*, 2019





# Principles of nascent RNA isolation

Table 2 | Methods used to investigate different steps of transcription


Method	Transcription step						
	TSS <sup>a</sup>	RNA capping	Promoter-proximal pausing	Co-transcriptional RNA processing	Transcription termination	Pol II CTD modification	Transcription bursting
<b>Chromatin isolation-based methods</b>							
caRNA-seq	No	No	No	Yes <sup>42,105–107</sup>	No	No	No
Start-seq	Yes <sup>43</sup>	No	Yes <sup>43</sup>	No	No	No	No
mNET-seq	No	No	Yes <sup>41,73</sup>	Yes <sup>41,63,64</sup>	Yes <sup>41</sup>	Yes <sup>41,63</sup>	No
SMIT-seq	No	No	No	Yes <sup>159,160</sup>	No	No	No
<b>Run-on methods</b>							
GRO-cap and PRO-cap	Yes <sup>4,42</sup>	No	No	No	No	No	No
GRO-seq, PRO-seq and ChRO-seq	No	No	Yes <sup>42,48,74</sup>	Yes <sup>166</sup>	Yes <sup>42</sup>	No	No
CoPRO	Yes <sup>49</sup>	Yes <sup>49</sup>	Yes <sup>49</sup>	No	No	No	No
<b>Metabolic labelling methods</b>							
TT-seq	No	No	No	No	Yes <sup>47</sup>	No	No
<b>Imaging-based methods</b>							
Intron sequential FISH	No	No	No	No	No	No	Yes <sup>55</sup>


caRNA-seq, chromatin-associated RNA sequencing; ChRO-seq, chromatin run-on and sequencing; CoPRO, coordinated precision run-on and sequencing; CTD, C-terminal domain; FISH, fluorescence in situ hybridization; GRO-cap, global run-on with cap selection; GRO-seq, global run-on sequencing; mNET-seq, mammalian native elongating transcript sequencing; Pol II, RNA polymerase II; PRO-cap, precision run-on with cap selection; PRO-seq, precision run-on sequencing; SMIT-seq, single-molecule intron tracking sequencing; TSS, transcription start site; TT-seq, transient transcriptome sequencing. "Yes" indicates that the method was designed specifically to detect TSSs, but GRO-seq can infer TSSs, as can PRO-seq and mNET-seq with sufficiently long reads.


# Intron seqFISH




















# Data analysis


- Spike-Ins
- Random nucleotides

# Why random nucleotides?



# Why random nucleotides?

- Adapter ligation



# Why random nucleotides?

- Adapter ligation



RNA with  
adapters



RT+PCR



# Data analysis work flow

- Eliminate duplicates (dedupe)



Diagram illustrating deduplication of data. It shows two groups of data, each consisting of two identical rows. Each row has a black bar, followed by 'NNNN' in blue, a red bar, and 'NNNN' in orange. The second group is identical to the first, with 'NNNN' in green instead of blue and orange.

Group 1 (Blue/Orange):

- Black bar, NNNN, Red bar, NNNN
- Black bar, NNNN, Red bar, NNNN

Group 2 (Green):

- Black bar, NNNN, Red bar, NNNN
- Black bar, NNNN, Red bar, NNNN



Diagram illustrating deduplication of data. It shows a single group of four identical rows. Each row has a black bar, followed by 'NNNN' in blue, a green bar, and 'NNNN' in black. The rows are stacked vertically, with each subsequent row shifted to the right.

Group 3 (Green):

- Black bar, NNNN, Green bar, NNNN
- Black bar, NNNN, Green bar, NNNN
- Black bar, NNNN, Green bar, NNNN
- Black bar, NNNN, Green bar, NNNN



# Data analysis work flow

- Eliminate duplicates (dedupe)
- Eliminate adaptors and random nmers (trim\_galore)
- Map reads to genome (bowtie2; merged genome from 2 species!)
- Separate reads by which genome and which strand they map to (samtools view -b and -f/-F 16; keep in mind that reads are reverse complement of actual transcripts)
- (optional: clip all but first base pair of every read, get Clip Reads)
- Normalise to Spike-In sequencing depth while transforming to .bedgraph file (genomeCoverageBed)
- Transform to .bw files (bedGraphToBigWig), visualise on genome browser
- Use featureCounts function in R to count reads in every gene and do e. g. differential expression analysis (manually divide resulting counts by Spike-In sequencing depth)

